

HPLC Column Storage Cabinet

SPECIFICATIONS:

- Capacity: 125 Columns
- Dimensions (WxHxD) : 280 X 1210 X 410mm
- No. of Trays: 5 Trays
- No. of Columns in a Tray: 25 Columns
- Stainless Steel 316 Quality
- Individual Lockable Tray
- Rust Proof



SPECIFICATIONS:

- Capacity: 500 Columns
- Dimensions (WxHxD) : 700 X 1400 X 700mm
- No. of Trays: 10 Trays
- No. of Columns in a Tray: 50 Columns
- Stainless Steel 316 Quality
- Individual Lockable Tray
- Rust Proof



Website

