

## HPLC Column Washing Pump

### SPECIFICATIONS:

- Flow Rate: 0.1 to 10.00 mL/min for  
10 mL/min head
- Pressure: 0 to 6,000 psi
- Flow Accuracy:  $\pm 5\%$  for a Flow rate of 0.1 to  
10.000 mL/min  
with 100% water @1000 psi
- Dimensions: 5.5" High x 10.375" Wide x 17.5" deep
- Weight: 24 lb
- Power: 100-120 VAC, 50-60 Hz; or  
220-240 VAC, 50-60 Hz
- Features: 4 Line LCD Display, Splash  
Proof Keypad  
Prime Purge Valve, Auto Stop  
Timer, Service Mode
- Optional: Pulse damper
- No. of Column: Upto 4 or 5 columns can be washed  
simultaneously
- Remote Inputs: RS-232



Website

